

**Camp Medica D.** - No. 29 Stanei Street, S4, Bucharest, Romania  
 phone: +4021-450 58 90  
 e-mail: [export@campmedica.ro](mailto:export@campmedica.ro)  
 http : //www.campmedica.ro



## “See Now” Coronavirus 2019 nCoV antibody Cassette Test

Whole blood/Serum/Plasma

For in vitro Diagnosis Use  
 Product Code: SN 8.5

### INTENDED USE

The “See Now” Coronavirus 2019 nCoV Ab. test is used for the qualitative determination of the novel coronavirus (SARS-CoV-2) antibody (IgG/IgM) in human serum, plasma or whole blood.

### PRINCIPLE

The “See Now” Coronavirus 2019 nCoV Ab. test is a chromatographic immunoassay (CIA) for the rapid qualitative determination of human novel coronavirus (SARS-CoV-2) IgG/IgM antibody in serum, plasma or whole blood. When the sample contains the SARS-CoV-2 antibody, it forms a complex with the gold label antigen SARS-CoV-2 (2019 nCoV antigen-colloidal gold conjugate). The complex moves forward under the action of chromatography and combines with the coated antibody (Mouse anti-human IgG and IgM monoclonal antibody) at the T lines (G and M) to form a complex and develop color ( lines G and M), which is a positive result. When the sample does not contain the SARS-CoV-2 antibody, no complex can be formed at the G and M lines and no red band appears, which is a negative result.

A build-in control line (C) always will appear in the test window when the test has performed properly regardless of the presence of absence of SARS-CoV-2 antibody in the specimen.

### MATERIALS SUPPLIED

- Test Device, capillary pipette, desiccant, Test Instruction, Sample buffer (PBS)

### MATERIALS REQUIRED BUT NOT SUPPLIED

- Specimen collection container, Timer

### STORAGE AND STABILITY

- The “See Now” 2019 nCoV IgG/IgM Ab Test should be stored at room temperature (4-30°C) in the sealed pouch or desiccated container.
- Do not use it after the expiration date.

### PRECAUTIONS

- FOR PROFESSIONAL AND IN VITRO DIAGNOSTIC USE ONLY.
- The test device should remain in the sealed pouch until use.
- There should be no smoking or eating where antigen containing materials are being handled. Wear disposable gloves and lab coat while handling specimens. Wash hands thoroughly afterwards.
- Decontaminate and dispose of specimens and all potentially contaminated materials if they contain infectious agent.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory have been evaluated.

### SPECIMEN COLLECTION AND STORAGE

- For serum samples collect blood in a tube without anticoagulant and allow it to clot.
- For whole blood or plasma samples collect blood in a tube containing anticoagulant (EDTA, citrate or heparin, respectively).
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
- **Finger puncture whole blood**

Clean the finger with an alcohol pad and let dry

Take a lancet and make a quick deep slab on the side of the finger.

Avoid squeezing the fingertip to accelerate bleeding as this tends to dilute the blood with excess tissue fluid.

- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature

### TEST PROCEDURE

- Remove the test device from pouch when ready to perform the test. Label the test device with patient or control identification
- Remove the test device from the sealed pouch by tearing at the notch. Then place the testing device on a level surface
- **For SERUM/PLASMA/ WHOLE BLOOD:** Holding the **capillary pipette** vertically, add **2µL** of specimen without air bubbles into the sample well, marked **S1**, then add **two (2) drops (80-100µl)** of sample buffer into the well marked **S**.
- Read the results in 15-20 minutes. Ensure that the background of the test area is white before interpreting the result. **Do not interpret result after 20 minutes.**

### INTERPRETATION OF RESULTS

#### Positive

Distinct pink colored bands (2 or 3) appear at the control and test lines regions (C, G, M lines).

#### Negative

Only one pink colored band appears at the control region (C).

#### Invalid

No visible band at the control region (C). Repeat with a new test device. If test still fails, please contact the distributor with the lot number.

### LIMITATION OF PROCEDURE

- The “See Now” 2019 nCoV Ab test device is used for the detection of SARS-CoV-2, **IgG/IgM Ab.** in human serum, plasma, whole blood. Based on a single reactive test result, a sample should not be considered **2019 nCoV IgG/IgM Ab.** positive. Further testing, including confirmatory testing, should be performed before a specimen is considered positive for **2019 nCoV IgG/IgM Ab.** Specimens containing precipitate may give inconsistent test results.
- The results of the reagent are only for clinical reference, which is not the only basis for clinical diagnosis and treatment. A confirmed diagnosis and treatment should only be made by a physician after all clinical and laboratory findings have been evaluated.

### QUALITY CONTROL

A procedural control is included in the test. A colored band appearing on the control region of the membrane indicates proper performance and reactive reagents.

### PERFORMANCE CHARACTERISTICS

- **Relative Sensitivity - 90%**
- **Relative Specificity - 99%**

