

**“See Now” CK-MB test  
 Serum/Plasma/Whole blood**



**For in vitro Diagnosis Use  
 Product Code: SN 6.3**

**Intended Use**

The “See Now” CK –MB test is an immunochromatography based one step in vitro test. It is designed for qualitative determination of CK-MB in human serum, plasma or whole blood specimens as an aid in the diagnosis of myocardial infarction.

**Principle**

Rapid CK-MB test is a sandwich immunoassay. When serum sample is added to sample pad, it moves through the conjugate pad and mobilizes gold anti-CK-MB conjugate that is coated on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-CK-MB antibody that is coated on the test region. If CK-MB is present at levels of 7.0 ng/ml or greater, the result is the formation of a colored band in the test region. If there is no CK-MB in the sample, the area will remain colorless. The sample continues to move to the control area and forms a pink to purple color, indicating the test is working and the result is valid.

**Precautions**

- **FOR IN VITRO DIAGNOSTIC USE ONLY**
- Do not use it after the expiration date.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.
- Handle all specimens as potentially infectious.
- Humidity and temperature can adversely affect results.

**Materials Provided**

- Test devices
- Pipette
- An instruction insert

**Materials Required But Not Provided**

Specimen collection container Centrifuge (for plasma only) Timer

**Storage and Stability**

The kit can be stored at room temperature (15-25°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.**

**Specimen Collection and Preparation**

- The “See Now” CK-MB Test can be performed using either serum, plasma or whole blood.
- Separate the serum or plasma from blood as soon as possible to avoid hemolysis. Only clear, non-hemolyzed specimens can be used. Remove the serum or plasma from the clot or red cells, respectively, as soon as possible to avoid hemolysis. Lipemic, icteric, or hemolyzed specimens may give inconsistent test result. Specimens, of containing precipitate, should be clarified prior to testing.


- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 4-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with national regulations covering the transportation of etiologic agents.

**Test Procedure**

Allow the test device, buffer, specimen to equilibrate to room temperature (15-25°C) prior to testing.

1. Remove the test device from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test device on a clean and level surface. Hold the dropper vertically and transfer 3-4 full drops (**120 - 160 µL**) of serum, plasma or whole blood. Avoid air bubbles.
3. Wait for the red line to appear. The result should be read at **10 minutes. Do not interpret the result after 30 minutes.**

**Interpretation of Results**



**Positive**

Distinct pink colored bands appear at the control and test line regions

**Negative**

Only one pink colored band appears at the control region.

**Invalid**

No visible band at the control region. Repeat with a new test device. If test still fails, please contact the distributor with the lot number.

**Quality Control**

A procedural control is included in the test. A red line appearing in the control region (C) is considered as an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

**Detection Limit**

“SeeNow” CK-MB, can detect Creatine Kinase with concentration of **7.0 ng/mL** or greater.